



Hillrom™



# TRUPOINT® SYSTEM SUPPORT HEAD

THINK LATERALLY

## IDEAL CONDITIONS FOR ANESTHESIA WORKSTATIONS WITH THE HORIZONTAL TRUPOINT SYSTEM SUPPORT HEAD

The reliable TruPort system modularity combined with space saving horizontal orientation makes TruPort System Horizontal the ideal solution for the anesthesia workstation. Positioned directly above the anesthesia machine, TruPort System Horizontal provides more freedom of movement in the OR and enables efficient workflow around the anesthesia workstation. Through modular architecture, supply modules for gases, power and data can be simply repositioned or even extended at any time according to your requirements – even post installation.




## MORE FLEXIBILITY FOR ANESTHESIA

- 2 lengths available (500 mm / 19.7" and 750 mm / 29.5")
- To be combined with all TruPort® arm systems
- Up to 4 side rails laterally mounted
- Accessories can be attached with just one click
- Maximum pay load of the head is 80 kg / 176 lbs
- Suitable TruPort system accessories available



## MORE SPACE IN THE OR

Compact architecture of the horizontal TruPort system support head. Positioned directly below the ceiling allows maximum head-room and freedom of movement.



## FLEXIBLE SOLUTIONS FOR THE FUTURE

React quickly and efficiently to new requirements in the OR. The simple repositioning and possibility to retrofit supply modules for gases, power and data, TruPort System Horizontal provides the flexible solution for the future.



## USER-FRIENDLY CONTROL

User-friendly control handles and an optional remote control allow ergonomic control for the user. Compatible with all arm systems, TruPort System Horizontal can be easily controlled in all situations.

For more information please contact your local Hillrom representative.

**hillrom.com**

Hillrom  
Amerika Building (7th Floor), Hoogoorddreef 15  
1101 BA Amsterdam, The Netherlands

TruPort naming corresponds to the following CE registered products according to MDD: TruPort 5000; TruPort 5500; TruPort 7000; TruPort 7500; Pre-Install Set - 4036100

The TruPort is intended to be used by clinicians and medically qualified personnel. Class IIb Notified Body: TÜV Rheinland LGA Products GmbH (CE0197), legal manufacturer TRUMPF Medizin Systeme GmbH + Co. KG, Carl-Zeiss-Straße 7-9, 07318 Saalfeld, Germany). This medical device is a regulated healthcare product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices. Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products. TRUMPF Medizin Systeme GmbH + Co. KG is a subsidiary of Hill-Rom Holdings, Inc.