



Welch Allyn
KleenSpec® Disposable Vaginal Specula
and Cordless Illumination System



An advancement in women's health
technology that's nothing short of brilliant

Welch Allyn KleenSpec® Disposable Vaginal Speculum and Cordless Illumination System—

now in every size

for every type of patient.



590 Series KleenSpec Disposable Vaginal Speculum

Enhance Patient Comfort

- The new extra small speculum allows you to perform a comfortable exam with virginal/pediatric patients as well as post-menopausal and post-hysterectomy women—no more referrals or uncomfortable exams using metal specula or the wrong size
- Smooth, molded acrylic designed to prevent pinching or binding of tissue for greater patient comfort

Help Reduce Cross-Contamination

- Single-use speculum helps to reduce cross-contamination for improved patient safety while eliminating the time and expense of reprocessing metal specs

Improve Staff Productivity

- Features a wider handle design for improved maneuverability and better balance during exams
- A strong grip and a yoke that's more flexible than low-cost substitutes



#59004
Large

#59001
Medium

#59000
Small

#590XS
Extra Small

Look for the distinctive Welch Allyn blue, green, white and yellow colored handles

“This light is brighter than any I have ever seen and shows true tissue color. I would recommend this system to any clinicians looking to improve the quality of their care.”

Shelley Kolton, MD / Downtown Women OB/Gyn Associates / New York City, NY



KleenSpec Cordless Illumination System

- LED technology offers a whiter, brighter light and true tissue color during an exam
- No bulb replacement necessary for a reduced cost of ownership
- Features a cordless design so you'll never have to worry about cords getting in your way, breaking or needing to be cleaned, speeding the exam and giving you a full range of motion
- Fits every size speculum
- Uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- Produces much less heat than traditional lamps and has a continuous on time of 80 minutes
- A lithium ion battery provides a recharge time of five hours



A site that conducts 75 exams per week can **save more than \$8,000 per year** by switching from metal to disposable specula*

*Based on estimated costs of specula, training, and salary. Estimated cost of metal vaginal specula from <http://gynexcorporation.com/stainlesssteel-2.aspx>. Salary amounts based on data from <http://www.bls.gov/oes/2011/may/oes291111.htm>

Ordering Information

KleenSpec Vaginal Specula—Premium 590 Series

- 590XS KleenSpec 590 Series Premium Disposable Vaginal Specula, **X-Small (24/box; 4 boxes/case)**
- 59000 KleenSpec 590 Series Premium Disposable Vaginal Specula, **Small (24/box; 4 boxes/case)**
- 59001 KleenSpec 590 Series Premium Disposable Vaginal Specula, **Medium (24/box; 4 boxes/case)**
- 59004 KleenSpec 590 Series Premium Disposable Vaginal Specula, **Large (18/box; 4 boxes/case)**

KleenSpec Vaginal Specula—Premium 590 Series with Smoke Tube

- 59005 KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, **Small (12/box; 4 boxes/case)**
- 59006 KleenSpec 590 Series Premium Disposable Vaginal Specula with Smoke Tube, **Medium (12/box; 4 boxes/case)**

KleenSpec Cordless Illumination System

- 80010 Complete KleenSpec Cordless Illuminator with Charging Station
- 80015 Same as above, for export purposes—specify country and voltage; conforms to IEC 60601-1
- 80000 KleenSpec Cordless Illuminator
- 74010 Charging Station for Cordless Illuminator
- 59010 Disposable Sheath for Cordless Illumination System **(100/box; 5 boxes/case)**

KleenSpec Storage Systems

- 58014 KleenSpec Vaginal Speculum Dispenser

Green Series™ Exam Light IV—Another LED innovation from Welch Allyn

- 48810 GS Exam Light IV with Table/Wall Mount
- 48950 Mobile Stand for GS Exam Light IV
- 52640 Disposable Sheath for Exam Light III and GS Exam Light IV **(25/box, 5 boxes/case)**

Explore the entire Welch Allyn women's health line, including the KleenSpec Series and Cordless Illumination System online at welchallyn.com/womenshealth or call us at **800.535.6663**.



Welch Allyn, Inc.
4341 State Street Road
Skaneateles Falls, NY 13153 USA
(p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

