



Hillrom™



2021 Product Catalog

HEALTHCARE FURNITURE

Distributed by Hillrom

Recliner

FEATURES AND BENEFITS:

- Infinite position recline
- Independent back and ottoman recline
- Robust flip-up arms with locking in all positions up/down
- Side retractable accessory hook
- Caregiver ergonomic push bar on back
- Front and rear lock casters with toe lock/release casters in front
- Compatible with Hillrom™ overbed tables and patient lifts
- Structure composed of injection molded plastic and cast aluminum with no exposed wood
- Molded foam cushions
- Solid surface arm caps available in multiple colors
- Available in non-porous fabrics with minimal sewing seams
- Optional central lock and steer

Dimensions	W: 33.75" D: 35" H: 44.25" SW: 23" SD: 20" SH: 19" Arm Height: 26.25" Back Seat Height: 28" Floor Clearance: 4.5"
Weight Capacity	500 lbs
Caster Size	4"
Range of full recline	74°
Part Number	Non-Powered: P9180 Powered: P9181
Overall Weight	Non-Powered: 114 lbs Powered: 125 lbs



Powered Recliner

- Independent power recliner and ottoman operation, with home return feature
- One integrated USB port
- Battery back-up

Bariatric Recliner

FEATURES AND BENEFITS:

- Infinite position recline up to 45°
- Independent back and ottoman recline
- Side retractable accessory hook
- Caregiver ergonomic push bar on back with release button
- Front and rear lock casters with toe lock/release
- Compatible with Hillrom™ overbed tables and patient lifts
- Structure composed of steel and hardwood with no exposed wood
- Molded foam cushions
- Solid surface arm caps available in multiple colors
- Available in non-porous fabrics with minimal sewing seams
- Optional removable arm (left side only)



Compatible with Hillrom™ overbed tables



Side retractable accessory hook

Dimensions	W: 36" D: 36.5" H: 45" SW: 29.75" SD: 20" SH: 20.5" Arm Height: 26.25" Back Seat Height: 30"
Weight Capacity	880 lbs
Caster Size	4"
Range of full recline	45°
Part Number	Non-Powered: P9190 Powered: P9191
Overall Weight	Non-Powered: 175 lbs Powered: 185 lbs

Powered Recliner

- Independent power recliner/ottoman
- One integrated USB port
- Battery back-up

20" Highback and Lowback

FEATURES AND BENEFITS:

- Ergonomic lumbar support back design for comfort
- Non-marring glides
- Wall saver back design
- Structure composed of steel and hardwood with no exposed wood
- Solid surface arm caps and powder coated steel frame available in multiple colors
- Available in non-porous fabrics with minimal sewing seams
- Optional flex back

Dimensions	Lowback: W: 23.5" D: 24.25" H: 33.75" SW: 19.5" SD: 18" SH: 18.5" Arm Height: 25.125"
	Highback: W: 23.5" D: 24.25" H: 40.75" SW: 19.5" SD: 18" SH: 18.5" Arm Height: 25.125"
Weight Capacity	300 lbs
Part Number	P9120
Overall Weight	Lowback: 28 lbs Highback: 30 lbs




Folding Chair with Wall Mount

FEATURES AND BENEFITS:

- Hinged powder coated steel frame available in multiple colors
- Wall mounted hanging bracket matching steel frame color
- Non-marring glides
- Wall saver back design
- Available in seamless, non-porous fabrics
- Mounting hardware for bracket included



Dimensions	Chair: W: 20" D: 19" H: 32.6875" SW: 17.75" SD: 17.75" SH: 17.75"
	Bracket: W:10" D: 2.15" H: 2.06"
Weight Capacity	300 lbs
Part Number	P9111
Overall Weight	Chair: 22 lbs Bracket: 1 lb



2000+ FABRIC OPTIONS

See page 18

30" and 40" Lowback

FEATURES AND BENEFITS:

- Ergonomic lumbar support back design for comfort
- Non-marring glides
- Wall saver back design
- Structure composed of steel and hardwood with no exposed wood
- Solid surface arm caps and powder coated steel frame available in multiple colors
- Available in non-porous fabrics with minimal sewing seams

30" Lowback

Dimensions	W: 33.5" D: 24.25" H: 33.75" SW: 29.5" SD: 18" SH: 18.5" Arm Height: 25.125"
Weight Capacity	400 lbs
Part Number	P9130
Overall Weight	37 lbs

40" Lowback

Dimensions	W: 43.5" D: 24.25" H: 33.75" SW: 39.5" SD: 18" SH: 18.5" Arm Height: 25.125"
Weight Capacity	500 lbs
Part Number	P9140
Overall Weight	47 lbs



Ergonomic lumbar support back design for comfort

Structure composed of steel and hardwood with no exposed wood

Sleeper Chair

FEATURES AND BENEFITS:


- 2 integrated USB ports
- Easy pull-out mechanism for chaise or sleep surface
- Interior cubby for personal storage while sleeping (24"W x 6"H x 4.625"D)
- Open bottom frame with 4 casters (2 front locking)
- Solid surface arm caps available in multiple colors
- Optional swivel, solid surface tray table (right side sitting)
- Optional LED reading light (right side sitting)

35" Sleeper Chair

Dimensions	W: 35" D: 37" H: 36" SW: 27" SD: 20" SH: 18" Arm Height: 24.5" Chaise: D: 64.125" Sleep Surface: D: 72"
Weight Capacity	350 lbs
Part Number	P9135
Overall Weight	170 lbs

45" Sleeper Chair

Dimensions	W: 45" D: 37" H: 36" SW: 37" SD: 20" SH: 18" Arm Height: 24.5" Chaise: D: 64.125" Sleep Surface: D: 72"
Weight Capacity	350 lbs
Part Number	P9145
Overall Weight	220 lbs



2000+ FABRIC OPTIONS
See page 18



Integrated USB ports



Optional LED reading light



Interior cubby for personal storage while sleeping

Sleeper Sofa

FEATURES AND BENEFITS:

- Available in 76" or 86" width
- 2 integrated USB ports
- Easy fold down sleep pad
- Interior cubby for personal storage while sleeping (24"W x 6"H x 4.625"D)
- Moisture barrier on sleep surface
- Open bottom frame with 5 aluminum legs and non-marring glides
- Solid surface arm caps available in multiple colors
- Optional 5 casters with 2 front locking
- Optional 2 drawers for storage of pillows & blankets (24" W x 6" H x 22" D)
- Optional swivel, solid surface tray table (right side sitting)
- Optional LED reading light (right side sitting)

76" Sleeper Sofa

Dimensions	W: 76" D: 36" H: 33" SW: 67.5" SD: 21.5" SH: 17" Arm Height: 23.75" Sleep Surface: W: 67.5" D: 30"
Weight Capacity	750 lbs
Part Number	P9176
Overall Weight	175 lbs

86" Sleeper Sofa

Dimensions	W: 86" D: 36" H: 33" SW: 77.5" SD: 21.5" SH: 17" Arm Height: 23.75" Sleep Surface: W: 77.5" D: 30"
Weight Capacity	750 lbs
Part Number	P9186
Overall Weight	182 lbs



Art of Care® OBT636

FEATURES AND BENEFITS:

- Variable height adjustment
- Dual sliding food service tray
- Tuning fork base design for fit and function with beds & furniture
- Steel column provides durability
- Braille on all touch points
- Automatic raising tabletop with easy lift handle
- Easy roll casters
- Thermofoil table surface
- Optional low base to fit under 1.8" surfaces

Dimensions	Overall Height: Standard Base (hi/hi) 48.9" (low/low) 31.4" Low Base: hi/hi) 46" (low/low) 29" Overall Width: 20.9" Food Service Tray Extended Width: 36.3" Overall Length: 41" Base Height: Standard Base: 2.5" Low Base: 1.8"
Weight Capacity	25 lbs
Part Number	P636
Overall Weight	Standard Base: 66 lbs Low Base: 66 lbs



Style 1
Flush top



Style 2
Anti-drip perimeter groove



Style 1
Flush with mirror/hand
grip cut out



Style 3
Food tray indentation with
mirror/hand grip cut out and
cupholder indentations



Style 4
Flush with no cut outs
or indentation

Art of Care® OBT635

FEATURES AND BENEFITS:

- Variable height adjustment
- Tuning fork base design for fit and function with beds & furniture
- Steel column provides durability
- Braille on all touch points
- Automatic raising tabletop with easy lift handle
- Easy roll casters
- Thermofoil table surface
- Optional low base to fit under 1.8" surfaces

Dimensions	Overall Height: Standard Base (hi/hi) 46.3" (low/low) 29.3" Low Base: hi/hi) 46" (low/ low) 29" Overall Width: 19.5" Overall Length: 39.5" Base Height: Standard Base: 2.5" Low Base: 1.8"
Weight Capacity	25 lbs
Part Number	P635
Overall Weight	Standard Base: 59 lbs Low Base: 59 lbs

Art of Care® OBT634

FEATURES AND BENEFITS:

- Variable height adjustment from 27.3" to 39.1"
- "C" base design for fit and function with beds & furniture
- Automatic raising tabletop with easy lift handle
- Easy roll casters
- Thermofoil table surface

Dimensions	Overall Height: (hi/hi) 39.1" (low/low) 27.3" Overall Width: 18.1" Overall Length: 34.8" Base Height: 3.6"
Weight Capacity	25 lbs
Part Number	P634
Overall Weight	31 lbs



Style 1
Flush top



Metropolitan®
P528 / BSCM

Dimensions	Metropolitan®: W-22.875", D-16.875", H-31.5" Aero®: W-20", D-17.5", H-30"
Weight Capacity	25 lbs
Part Number	P528-Metropolitan P526-Aero
Overall Weight	92 lbs-Metropolitan 68 lbs-Aero

Art of Care® Bedside Cabinets

FEATURES AND BENEFITS:

- Hardwood and hardwood veneer drawer fronts and side panels for style
- Ventilated back
- Steel frame and base
- Easy roll casters
- Three drawer with easy grip drawer handles
- Ultra tough top, stain and scratch resistant
- Rounded drawer corners for ease of cleaning
- Available in multiple thermofoils and laminates



Aero®
P526 / BSCR

Fabric & Finish Studio

DESIGN YOUR OWN:


This feature gives you the ability to view fabric and finish offerings and "virtually" apply them to the product for easy viewing. You can also order fabrics and finishes at no charge.



FABRIC OPTIONS:

With nearly 2000 options, Hillrom's fabric options are as versatile as the needs of today's healthcare environments. Choose from a wide selection of coated and performance fabrics, ranging from faux leathers to biophilic designs.

HILLROM OFFERS MULTIPLE SOLID COLOR OPTIONS, THERMOFOIL, WOOD, AND POWDER COAT FINISH OPTIONS TO MATCH YOUR HOSPITAL'S DESIGN AND AESTHETICS.



Visit hillrom.com to "Design Your Own" furniture.

2000+
Fabric Options





Hillrom™

ABOUT HILLROM

Hillrom is a global medical technology leader whose 10,000 employees have a single purpose: enhancing outcomes for patients and their caregivers by advancing connected care. Around the world, our innovations touch over 7 million patients each day. They help enable earlier diagnosis and treatment, optimize surgical efficiency and accelerate patient recovery while simplifying clinical communication and shifting care closer to home. We make these outcomes possible through connected smart beds, patient lifts, patient assessment and monitoring technologies, caregiver collaboration tools, respiratory care devices, advanced operating room equipment and more, delivering actionable, real-time insights at the point of care. Learn more at hillrom.com.

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

130 E. Randolph St. Suite 1000, Chicago, IL 60601

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. APR176801 rev 3 15-JUL-2021 ENG – US