










DISPOSABLE CUFF SIZE, CONFIGURATION, PART NUMBER REFERENCE GUIDE

		NO TUBES OR CONNECTERS REQUIRED	One-Tube				Two-Tube			
										
Size (20 cuffs/case)	Range	FlexiPort "Raw Cuff"	Screw	Tri-purpose	Locking	Bayonet	Screw	Locking	Tri-Purpose & Empty Tube	Inflation Sys.
Small Infant (6)	7-10 cm	SOFT-06	SOFT-06-ISC	SOFT-06-1TP	SOFT-06-1MQ	SOFT-06-1HP	N/A	N/A	N/A	N/A
Infant (7)	9-13 cm	SOFT-07	SOFT-07-ISC	SOFT-07-1TP	SOFT-07-1MQ	SOFT-07-1HP	SOFT-07-2SC	SOFT-07-2MQ	SOFT-07-2TP	SOFT-07-2BV
Sm. Child (8)	12-16 cm	SOFT-08	SOFT-08-ISC	SOFT-08-1TP	SOFT-08-1MQ	SOFT-08-1HP	SOFT-08-2SC	SOFT-08-2MQ	SOFT-08-2TP	SOFT-08-2BV
Child (9)	15-21 cm	SOFT-09	SOFT-09-ISC	SOFT-09-1TP	SOFT-09-1MQ	SOFT-09-1HP	SOFT-09-2SC	SOFT-09-2MQ	SOFT-09-2TP	SOFT-09-2BV
Sm. Adult (10)	20-26 cm	SOFT-10	SOFT-10-ISC	SOFT-10-1TP	SOFT-10-1MQ	SOFT-10-1HP	SOFT-10-2SC	SOFT-10-2MQ	SOFT-10-2TP	SOFT-10-2BV
Adult (11)	25-34 cm	SOFT-11	SOFT-11-ISC	SOFT-11-1TP	SOFT-11-1MQ	SOFT-11-1HP	SOFT-11-2SC	SOFT-11-2MQ	SOFT-11-2TP	SOFT-11-2BV
Adult Long (11L)	25-34 cm	SOFT-11L	SOFT-11L-ISC	SOFT-11L-1TP	SOFT-11L-1MQ	SOFT-11L-1HP	SOFT-11L-2SC	SOFT-11L-2MQ	SOFT-11L-2TP	SOFT-11L-2BV
Lg. Adult (12)	32-43 cm	SOFT-12	SOFT-12-ISC	SOFT-12-1TP	SOFT-12-1MQ	SOFT-12-1HP	SOFT-12-2SC	SOFT-12-2MQ	SOFT-12-2TP	SOFT-12-2BV
Lg. Adult Long (12L)	32-43 cm	SOFT-12L	SOFT-12L-ISC	SOFT-12L-1TP	SOFT-12L-1MQ	SOFT-12L-1HP	SOFT-12L-2SC	SOFT-12L-2MQ	SOFT-12L-2TP	SOFT-12L-2BV
Thigh (13)	40-55 cm	SOFT-13	SOFT-13-ISC	SOFT-13-1TP	SOFT-13-1MQ	SOFT-13-1HP	SOFT-13-2SC	SOFT-13-2MQ	SOFT-13-2TP	SOFT-13-2BV

For vinyl cuffs, replace "SOFT" with "VINYL" in part number (Vinyl "TP" version not available).

Neonatal Soft Disposable Cuffs

Neonatal Monitors					
connector		SnapQuik™		Luer SLip	
Size (10 cuffs/box)	Range	One-Tube	Two-Tube	One-Tube	Two-Tube
Neo 1	3.3 - 5.6 cm	Neo-1-1	Neo-1-2	5082-101-1	5082-101-2
Neo 2	4.2 - 7.1 cm	Neo-2-1	Neo-2-2	5082-102-1	5082-102-2
Neo 3	5.4 - 9.1 cm	Neo-3-1	Neo-3-2	5082-103-1	5082-103-2
Neo 4	6.9 - 11.7 cm	Neo-4-1	Neo-4-2	5082-104-1	5082-104-2
Neo 5	8.9 - 15.0 cm	Neo-5-1	Neo-5-2	5082-105-1	5082-105-2

SnapQuik™ is a registered Trademark of Colder Products Company