



Hillrom™

PROTECTING PATIENT INFORMATION AND YOUR NETWORK

The ELI 380 ECG and Q-Stress System now have a Conditional Authority to Operate (ATO-C) on Department of Defense networks.

DOD RMF



APPROVED



Welch Allyn®
ELI® 380 ECG



Welch Allyn
Q-Stress®
Cardiac
Stress Testing
System

Connecting to an EMR or ECG management system promotes a streamlined, error-free workflow. When doing so, it's critical to consider the security of both your connected systems and the environments they connect to.

The Welch Allyn ELI 380 ECG and Q-Stress System include features that meet strict security requirements from the United States Department of Defense (DoD), so you can take advantage of connected workflows that help improve patient care and information access.

Simple. Secure. Connected.

Help protect patient privacy and data integrity while safeguarding your network. Enhanced security features help keep your data safe where it is and when it moves, so you can protect the information that matters most. Our secure solutions offer:

- FIPS 140-2 certified wireless radio*
- Encryption of data stored on the device or transmitted over the network
- Connect to the data management system of your choice through standards-based protocols like HL7® or DICOM®
- Controlled user access through LDAP or Active Directory—allowing for unique usernames and passwords, along with role-based access to PHI and settings

ELI® 380 ECG Workflow



Q-Stress® System Workflow



Contact your Welch Allyn representative or visit www.hillrom.com to learn more.

hillrom.com

4341 State Street Road, Skaneateles Falls, NY 13153

* ELI 380 ECG only.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2020 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR52702 Rev 1 20-MAY-2020 ENG – US