



Hillrom™

PERIODIC INSPECTIONS



Mobile, Freestanding & Overhead Lift Systems

REDUCE SAFETY RISKS AND STAY COMPLIANT WITH:

✓	Documentation <ul style="list-style-type: none"> Joint commission requirements
✓	Operational Inspection <ul style="list-style-type: none"> Sling bar Lift Strap Decal Emergency stop Wheels Lift arms Hand controls Brakes
✓	Load Testing <ul style="list-style-type: none"> Maximum load testing on all lift components Emergency lowering
✓	Environmental Impact <ul style="list-style-type: none"> Corrosive environments Safety critical fixing point inspection
✓	Additional Components for Overhead Rail System <ul style="list-style-type: none"> Rail switches Ultra twist Transfer motor traverse Activity rail Pool side lifter Check stops

For more information, please contact your Hillrom sales representative at 1-800-445-3720.

hillrom.com

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.