



Nearly 100 years in the making...

Welch Allyn 767 Wall and Mobile Aneroids

Improved aesthetics and a new modern design make the 767 Aneroid an easy-to-use and attractive blood pressure aneroid. Based on Welch Allyn's near century of experience building and designing innovative blood pressure technology, it's the smart choice for accurate and reliable readings in your facility.

WelchAllyn[®]

Advancing Frontline Care[™]



Benefits that matter.

The 767 Aneroid is a testament to Welch Allyn's nearly 100-year history of blood pressure innovation—combining the traditional benefits you have come to expect from Welch Allyn with improved aesthetics and a new, contemporary design.

Accurate

- Laser-engraved dial face provides unmatched accuracy
- Certified accurate to ± 3 mm Hg
- Recessed dial face and high-contrast red pointer increase visibility and reduce parallax errors

Reliable

- Jeweled movement contributes to a longer life
- Lifetime calibration warranty proves our commitment to durability

Convenient

- New 5-leg stand provides stability and maneuverability
- Built-in storage basket on wall and desktop model

Safe

- FlexiPort® blood pressure cuffs are antimicrobial treated to help prevent mildew and bacteria growth
- Latex-free
- Mercury-free

Ordering Information

767 Wall Aneroid

7670-01	767 Wall Aneroid with Reusable One-Piece Adult Cuff*
7670-01C	767 Wall Aneroid with Reusable One-Piece Child Cuff
7670-01KPA	767 Wall Aneroid (KPA) with Adult Cuff
7670-02	767 Wall Aneroid with 8 ft. Coiled Tubing (No cuff included)

767 Mobile Aneroid

7670-03	767 Mobile Aneroid on 5-Leg Stand with Reusable One-Piece Adult Cuff*
7670-03K	767 Mobile Aneroid on 5-Leg Stand Only (No cuff included)
7670-10	767 Mobile Aneroid with Reusable One-Piece Adult Cuff*
7670-10K	767 Mobile Aneroid with Reusable One-Piece Child Cuff
7670-10KPA	767 Mobile Aneroid (KPA) on 5-Leg Stand with Reusable One-Piece Adult Cuff
7670-70	767 Aneroid with Mobile Bracket

767 Desktop Aneroid

7670-16	767 Desktop Model with Reusable One-Piece Adult Cuff
7670-16KPA	767 Desktop Model (KPA) with Adult Cuff

*Model available with traditional two-piece cuff and bladder inflations system.
Add "CB" to the end of part number (e.g. 7670-01CB).



Simple color changes and new rounded edges give the aneroid a modern look.



A new 5-leg mobile stand provides stability and maneuverability.

Not all models are available in all markets, check with your local Welch Allyn Representative. For more information about our blood pressure solutions, visit welchallyn.com.

WelchAllyn®

Advancing Frontline Care™

Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220
Skaneateles Falls, NY 13153-0220 U.S.A.
Telephone: +1.315.685.4602
+1.800.535.6663
Fax: +1.315.685.0054
E-mail: info@welchallyn.com

Welch Allyn European Customer Service: +353 46 90 67790
Australia: +61 29 638 3000
China: +86 21 6327 9631
France: +33 155 69 58 49
Germany: +49 6950 985 132
Latin America (Miami office): +1 305 669 9003

Malaysia: +603 78753341
Netherlands: +31 202 061 360
Singapore: +65 6419 8100
South Africa: +27 11 777 7555
Sweden: +46 858 536 551
United Kingdom: +44 207 365 6780