



Welch Allyn CP 50™ Electrocardiograph



Simple. Affordable. Compact.

The practical CP 50 ECG—designed to improve office efficiencies.


- > Offered at an affordable price to help maximize your budget
- > Easy-to-use, with one-button operation, makes performing ECG tests efficient
- > High-resolution thermal printer provides quick outputs
- > Full-color, touch-screen display to help improve staff satisfaction
- > Battery operation lets you easily transport your device to the patient
- > Optional Mobile Stand and Carrying Case are also available
- > Optional adult and pediatric MEANS interpretation offers a second opinion



EASY TO USE

Touch-Screen Display


Facilitates one-button operation, guiding you through each test quickly and easily.



EASY TO FOLLOW

On-Screen Instructions

Saves valuable time and helps ensure accurate data capture by highlighting poor signals and connections.



EASY TO READ

Programmable Report Formats

Prints reports in the format that works best for you with a high-resolution thermal printer.

WelchAllyn®

Advancing Frontline Care™

Improve your practice workflow and efficiency with the Welch Allyn CP 50™ ECG



Lightweight and portable

- One-button operation facilitates quick setup and operation
- Optional carrying case
- Optional mobile stand
- Weighs just over 4.4 lbs



Welch Allyn CardioPerfect™ Workstation

Bring your cardiopulmonary data together with the Welch Allyn CardioPerfect® Workstation

- Manage patient information and tests performed with the full suite of Welch Allyn cardiopulmonary devices
- Get paperless efficiency
- Store, access, review, compare and edit diagnostic test data from any computer on your network
- Retrieve records with the touch of a key, and review and edit right on the screen
- Access information from your exam rooms, your office or at home



CP 50 Rolling Cart, Part # 406815

Ordering/Parts information for the CP 50

Device configurations:

- CP50-1EN1** CP 50—Non-interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.
CP50A-1EN1 CP 50—Interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.
CP50P-1EN1 CP 50 Plus—Non-interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.
CP50AP-1EN1 CP 50 Plus—Interpretive 12-lead multichannel ECG, color touch-screen display and 4" thermal Printer.

Accessories:

- | | | | |
|-------------------|---|---------------|--|
| 715006 | ECG Multi-function Electrode Adapter | 406020 | CP 50 Printer Paper, Roll 4 rolls/case |
| 45008-0000 | Resting Tab Electrodes (box of 1000) | 715107 | CP 50 Carrying Case |
| 719653 | 10-lead Patient Cable, AHA, Banana, CP 50 | 406815 | CP 50 Cart |
| 406021 | CP 50 Printer Paper Z fold 4 packs/case | | |

Partners in CareSM Service Options:

- S1-CP50** CP 50, Comprehensive Partner Program
S1-CP50-2 CP 50, Comprehensive Partner Program 2YR

Warranty information

- CP 50 device—1 year Patient Cable /Battery—90 days

Flexible configuration options to meet the needs of your practice

Choose the CP 50 that's best for you



Welch Allyn CP 50™

- Reliable touch-screen display
- Optional MEANS interpretation
- High-resolution thermal printer
- Internal storage

Welch Allyn CP 50™ Plus with Connectivity

- Includes all of the features of the CP 50, plus:
- Flexible connectivity options
 - External printing capability
 - Programmable patient entry fields
 - Administrative controls
 - Store ECGs to a USB memory stick as PDF's

WelchAllyn®

Advancing Frontline Care™

Contact your local Welch Allyn or distributor representative for a full list of configurations, or visit www.welchallyn.com/CP50 today for more information.

Welch Allyn Corporate Headquarters

4341 State Street Road, Skaneateles Falls, NY 13153-0220 USA

Telephone: +1.800.535.6663

www.welchallyn.com/CP50

© 2013 Welch Allyn SM4098 Rev A