

WelchAllyn



Welch Allyn Connex® Vital Signs Monitor 6000 Series

Experience accurate, automatic spot-check vitals and continuous monitoring sent to your EMR from the bedside to help support improved patient safety and clinical decisions.

The Connex Vital Signs Monitor represents the next generation in vital signs monitoring with an easy-to-use color touchscreen and advanced vital signs technology to monitor vitals, respiration and motion. Plus, EMR integration gives clinicians immediate access to the information they need—where and when they need it.

Connect to the future of vital signs today

Designed specifically for busy floors like yours to increase efficiency

- > Captures both spot-check vitals and continuous monitoring
- > EMR integration for accurate electronic vitals documentation
- > Vivid color touchscreen simplifies training and operation
- > Comprehensive point-of-care documentation including customizable data fields
- > Advanced Patient ID including two-way name confirmation
- > Patient data protection for HIPAA compliance
 - Protect patient identifiable information with customizable options like abbreviations
 - Allow automated deletion of patient data after sending to EMR
 - Apply rules to check record completeness before saving
- > Easy-to-understand error messages for simple troubleshooting

Industry-leading parameters for accurate patient information

- > SureBP® blood pressure technology
 - Adult, pediatric and neonatal blood pressure modes
 - Customizable blood pressure monitoring programs help nurses follow hospital protocol
 - 15-second readings for greater patient comfort
- > SureTemp® Plus and/or Braun ThermoScan® PRO 6000 ear thermometry (optional—choose one or both)
- > Masimo® or Nellcor® SpO₂ algorithms (optional)—both sensors and signal processing
- > Masimo noninvasive continuous monitoring of total hemoglobin (optional)
- > Manual measurements entry (weight, height, respiration rate, pain) and documentation of modifiers like body position and O₂ therapy
- > Customizable documentation fields to fully document at the point of care




Upgradeable design to help protect your investment

The Connex Vital Signs Monitor is designed for the med/surg environment with an easy-to-use, configurable touchscreen display. And you can future-proof your purchase with these upgrade options designed to help keep patients safe and improve outcomes:

- > Braun ThermoScan PRO 6000 Ear Thermometer
- > Health o meter® Digital Weight Scale connectivity
- > Covidien® Microstream® Capnography (etCO₂)
- > Masimo® Acoustic Respiration (RRa)
- > Proven EMR connectivity
- > EarlySense® contact-free monitoring
- > Masimo® SpHb noninvasive hemoglobin
- > Best-in-class wireless technology


Advanced accessory management to simplify your life

- > Optimal organization of accessories
- > Helps prevent tangling and sensor damage
- > Cleaning supply storage encourages proper hygiene
- > LED lighting helps in dark rooms without disturbing patients



Protect your investment with service programs and tools

- > Flexible Partners in CareSM programs cover all your service needs
- > Fully customizable, including programmable location identifier and quick USB configuration cloning
- > PartnerConnectTM-ready for software-driven repairs and upgrades



Configurations Available

CONNEX VSM	6700	6800
Parameters		
SureBP NIBP	●	●
SureTemp Plus Thermometry	●	●
Braun PRO 6000 Ear Thermometry	●	●
Masimo or Nellcor SpO ₂	●	●
EarlySense Contact-Free Motion	↑	↑
Covidien Microstream EtCO ₂	↑	↑
Masimo RRa	↑	↑
Masimo SpHb	↑	↑
Communications		
4 USB Ports	●	●
Ethernet	●	●
802.11 a/b/g Wireless Radio	↑	●
Accessories		
Barcode Scanner	●	●
Weight Scale	●	●
EMR Connectivity	●	●
Printer	●	●

● Included
 ● Optional
 ↑ Upgradeable

ACCESSORIES

4800-60	Accessory Cable Management mobile stand with basket
4701-62	Wall mount with basket
4900-60	Accessory Cable Management Stand with basket, extended housing
008-0834-01	GCX wall mount
4900-62	GCX wall mount, extended housing
06000-900	Braun PRO 6000 dock with thermometer
6000-NC	Nurse call cable
6000-916HS	Medical 2D barcode scanner with coiled USB
6000-920	Wireless radio kit—internal 802.11 a/b/g
4500-925	USB cable for wired connectivity
660-0138-00	Ethernet patch cable 5 ft.
6000-BT0	One-year service contract, NIBP and temperature
6000-BTS	One-year service contract, NIBP, temperature, SpO ₂
S1-6000	Connex VSM 6000 Comprehensive Partnership Program, 1-year agreement
S1-6000-2	Connex VSM 6000 Comprehensive Partnership Program, 2-year agreement
S2-6000	Connex VSM 6000 Biomed Partnership Program, 1-year agreement
S2-6000-2	Connex VSM 6000 Biomed Partnership Program, 2-year agreement



For more information about the Connex Vital Signs Monitor, please contact your local Welch Allyn representative, call 1.800.535.6663 to speak with a customer service representative, or visit www.welchallyn.com.

Welch Allyn Corporate Headquarters
 4341 State Street Road, P.O. Box 220, Skaneateles Falls, NY 13153-0220 USA
 (p) 800.535.6663 (f) 315.685.3361

WWW.WELCHALLYN.COM

© 2017 Welch Allyn MC14317 80021073 Ver A 2017-07

