### **I·I'** Hillrom™



# SIMPLE. SECURE. CONNECTED.

The Welch Allyn® Q-Stress® Cardiac Stress Testing System brings proven clinical excellence to the forefront of diagnostic cardiac stress testing.

Experience the simplicity of our intuitive user interface and best-in-class algorithms so you can quickly access diagnostic information and exceptional accuracy in real-time ST-segment monitoring, arrhythmia detection and resting ECG interpretation.





The Welch Allyn Wireless Acquisition Module (WAM™) makes exams easier and creates a more comfortable, trip-free experience for your patients. Transfer patients to a Stress Echo table with ease and free your care area from unnecessary wires with a simple, untethered workflow that helps minimize artifact while capturing quality exam data.

#### **SIMPLE**

Left-to-right navigation enables quick access to the diagnostic information needed to make the best clinical decisions

Source Consistency Filter (SCF) helps minimize noise and baseline artifact

Real-time full disclosure allows clinicians to examine and label historical ECG data at any point in the exam

Best-in-class VERITAS®
algorithm offers exceptional
accuracy in real-time STSegment Monitoring, arrythmia
detection, and resting ECG
interpretation

The Q-Stress System is a safe, easy-to-use system. The large monitor provides real-time results that have proven to be more accurate than other systems we have used.

#### **Tiffany Piro**

CMRP Procurement Specialist Healthcare Associates of Texas

#### **SECURE**

Encryption helps protect the information that matters most

Security features help you take advantage of connected workflows for improved patient care and information access

Flexible user configurations tailor over-reads, billing, and data storage to your workflow

#### **CONNECTED**

Send data to virtually any EMR, ECG management system or PACS with our customizable, scalable platform

Communicate bidirectionally using HL7® or DICOM® to reduce errors while providing an electronic, editable format to your team



## CHALLENGE ACCEPTED

The Q-Stress® Solution can help address the clinical challenges you face with a simple workflow and clear data presentation.

Our complete clinical data presentation allows you to see patient information on one screen with a 24-inch touch or non-touch display. An intuitive left-to-right workflow leads you through stress tests with all settings just one click away.

#### 1. VITALS DASHBOARD

Visualize key vital signs—like HR, target HR, BP, SpO2, ST level and slope—to assess a patient's status

#### 2. CONTINUOUS QRS AVERAGES

Track ST changes, with current and reference data, through the superimposed 12 averaged ECG complexes

#### 3. FULL DISCLOSURE


Provides a single-lead view, highlighting events and isolating them for further evaluation complexes

#### **8. ACQUIRE AND DOCUMENT**

Access common functions at the touch of a button for 12-lead ECG, Averages, Rhythm Print, etc

#### 7. SIMPLE NAVIGATION

Intuitive menus lead you through test stages in an organized, left-to-right workflow



#### 4. FLEXIBLE ECG DISPLAY

Customize your view with one click, including displayed lead layout, gain, speed and filter settings

#### 5. ST DISPLAY

View graphical trends including ST changes in the 12 leads, plus a superimposed, expanded averaged complex

#### **6. GRAPHIC TRENDS DISPLAY**

Visualize live HR, Mets, BP and ST changes

# MAKING THE COMPLEX SIMPLE

We're helping you streamline workflows and promote data security with simple EMR connectivity. Install Q-Stress® software on a server and connect to multiple device workstations throughout your network to schedule exams, review full disclosure exams, enter conclusions and electronically sign results.



## SECURE YOUR DATA AT REST AND IN MOTION

- Help make security simple with features that meet strict requirements from the U.S. Department of Defense's Risk Management Framework (DoD RMF).
- The Q-Stress System supports your ability to save a PDF file to a common folder while utilizing networked workstations and review stations to share and access patient data securely.

## A SECURE, CONNECTED WORKFLOW ALLOWS YOU TO:

- Reduce costs by eliminating ECG paper
- Retrieve a stress test worklist that includes patient demographics and order information
- Minimize manual steps and transcription errors
- Access and send stress reports within your Healthcare IT System



**EHR OR PACS** 

Stress Test Worklist
DICOM or HL7



Q-STRESS®

Stress Test Report
DICOM, HL7 or XML with PDF



**EHR OR PACS** 

### CONNECT MEANINGFULLY

Electronic scheduling capabilities allow you to monitor and manage workflow tasks continuously and effectively. Take advantage of a true, secure paperless workflow where final reports, full disclosure data and all ECG events can be reviewed and edited electronically.

**APPROVED** 





#### **ORDERING INFORMATION**

| Part Number | Description                                                                       |  |
|-------------|-----------------------------------------------------------------------------------|--|
| QS6-ATTD1   | Pro Level Q-Stress® System with Touch Monitor and Treadmill; DICOM Connectivity   |  |
| QS6-ATTDX   | Pro Level Q-Stress System with Touch Monitor and no Treadmill; DICOM Connectivity |  |
| QS6-MLTC1   | Mid Level Q-Stress System with LCD Screen and Treadmill; PDF/XML Connectivity     |  |
| QS6-MLTCX   | Mid Level Q-Stress System with LCD Screen and no Treadmill; PDF/XML Connectivity  |  |
| 4100-036-51 | WAM™ Kit with AHA Short Lead Set                                                  |  |
| 9100-026-60 | Thermal Paper, 8.5 in x 11 in                                                     |  |

#### **SMARTCARE SERVICE PROGRAM OPTIONS**

| Program     | Product                     | Part Number      | Program Terms |
|-------------|-----------------------------|------------------|---------------|
| Protection+ | Q-Stress System w/Treadmill | S9-QSTM-PROPL-PS | 3 Years       |
|             | Q-Stress System             | S9-QS-PROPLUS-PS | 3 Years       |

#### **DELIVERY ASSIST**

Get up and running quickly with complimentary dock-to-door delivery, unboxing and on-site system preparation. Included with every Q-Stress purchase at no cost, Delivery Assist gets you up and running quickly with complimentary dock-to-door delivery, unboxing and on-site system preparation.

Experience simple, secure, connected stress exams. Learn more at hillrom.com/qstress today.

#### hillrom.com

4341 State Street Road, Skaneateles Falls, NY 13153